Prior to use, carefully read the instructions for use. **EN**

G-CEM LinkAce®

SELF-ADHESIVE RESIN CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges.
- Cementation of metal, ceramic, fiber posts, and cast post and cores.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

A. CEMENTATION TECHNIQUE for ceramics, indirect composites and metal-based restorations

- TOOTH PREPARATION**
 - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
 - Clean the prepared tooth with pumice and water. A clean surface is paramount for an optimal adhesion.
 - Rinse thoroughly with water. Dry thoroughly by gently blowing with oil free air. Prepared surface should be dry.

NOTE:

 - Do not use substances like H₂O₂, EDTA and/or disinfectants as these may lower the bond strength.
 - Special care should also be taken to remove so called imaging powders very thoroughly with pumice and water.
- RESTORATION PREPARATION**

Make sure that the restoration is pretreated and handled according to the manufacturer’s instructions.
- DISPENSING**
 - Before dispensing the first time from a new syringe or after a long interval in between use, bleed the pastes from the syringe.
 - Remove the cap of the syringe by rotating 1/4 turn anti-clockwise. Attach a GC Automix Tip by aligning the V shaped notch on the rim of the mixing tip with the V shaped notch between the syringe barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of the mixing tip 1/4 turn clockwise.
 - Extrude material directly into the restoration.

NOTE:

The used mixing tip can be left on the syringe to serve as a storage cap until next use. For each application, place a new mixing tip.
- CEMENTATION**
 - Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 45 seconds from start of mix at 23°C (73.4°F). Higher temperatures will shorten the working time.

NOTE:

Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.
 - Maintain moderate pressure until final polishing is done.
- EXCESS CEMENT REMOVAL**
 - Maintain moderate pressure and make sure the restoration remains in place.
 - To remove excess cement, light cure using a visible light curing unit for 1 second or allow the cement to self-cure until the excess cement feels rubbery.
 - Remove excess cement with appropriate instruments, holding the restoration in position.

- FINAL SET**

While maintaining moderate pressure, light cure all surfaces / margins for 20 seconds each (Halogen/LED, intensity: 700mW/cm², covering wavelength: 430 to 480 nm).

In case of the following restorations, let the material set for 4 minutes after removing excess cement.

 - Metal and metal coping restorations
 - Metal-free dark and opaque restorations
 - Translucent metal-free restorations greater than 2mm in thickness
- FINAL POLISHING**

If needed, finish and polish margins with appropriate instruments.

B. CEMENTATION TECHNIQUE for metal, ceramic, fiber posts, and cast post and cores

- TOOTH PREPARATION**

Prepare the post space according to manufacturer’s instructions. The post space should be cleaned, rinsed and dried thoroughly using paper points. The post space should be dry.

NOTE:

 - Sodium hypochlorite solution is recommended to use for the chemical cleaning of the post space.
 - Do not use H₂O₂ and/or EDTA to chemically clean the post space as these products may lower the bond strength.
- POST PREPARATION**

Make sure that the post is pretreated and handled according to the manufacturer’s instructions.
- DISPENSING & CEMENTATION**
 - Before dispensing the first time from a new syringe or after a long interval in between use, bleed the pastes from the syringe.
 - Attach a GC Automix Tip for endo. For details, see A 3b) above.
 - Apply the extension tip to the mixing tip.
 - Extrude material directly into the post space.
 - Insert the post immediately into the post space within 1 minute after cement application.

NOTE:

 - The oral temperature accelerates the setting reaction of the cement.
 - Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.
 - Keep the tip immersed in the cement while dispensing to avoid trapping air.
 - DO NOT use lentulo spinals to place the cement in the post space as this may shorten the working time.
 - Maintain moderate pressure.
- EXCESS CEMENT REMOVAL**

Continue to maintain moderate pressure making sure the post remains in place, and remove excess cement with appropriate instruments. In case of a cast post and core, tack cure all surfaces for 1 second each with a visible light curing unit.
- FINAL SET**
 - While maintaining moderate pressure, light cure all surfaces/margins for 20 seconds each (Halogen/LED, intensity: 700mW/cm², covering wavelength: 430 to 480 nm).
 - Let the material set for another 4 minutes.

STORAGE

Store at room temperature away from heat, moisture and direct sunlight (4-25°C) (39.2-77.0°F). (Shelf life: 2 years from date of manufacture)

SHADE

A2 (Vita® shade), A03 (opaque), B01 (opaque) and Translucent Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGE

- G-CEM LinkAce syringe 4.6 g (2.7 mL) (2), GC Automix Tip Regular (15), GC Automix Tip for endo with extension tip (5)
- GC Automix Tip Regular (20)
- GC Automix Tip Regular (10)
- GC Automix Tip for endo with extension tip (10)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Do not use G-CEM in combination with eugenol containing materials as eugenol may hinder G-CEM from setting or bonding properly.
- Use a protective light shield or similar protective eye wear during light curing.
- This product is not indicated for filling or core build-up.

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